

## CONVERSION TABLE FOR WALL MOUNTED OUTLETS

	OHIO MEDICAL/OHMEDA		PURITAN BENNETT			NCG/CHEMETRON		OXEQIP/MEDSTAR	
	<a href="#">DIAMOND 1</a>	<a href="#">DIAMOND 2 &amp; 3</a>	<a href="#">SERIES 73</a>	SERIES "A"	LITHARGE	SERIES 148,400,460	<a href="#">SERIES 500, 560</a>	SERIES # 7 OR 5/8-24	<a href="#">SMOOTH BORE OR 5/8-27</a>
INSTALLATION TOOL	O-VC-OH1-TOOL	Not needed	O-CPBTOOLKIT	Not needed	Not needed	O-VC-OH1-TOOL	Not needed	O-VC-OH1-TOOL	Not needed
KIT NUMBERS		<b>Complete Conversion</b>							
QD	O-OH1-QD-L-G	O-OH2-QD-L-G	O-C73-QD-L-G-W	O-CSA-QD-L-G-W	O-CLI-QD-L-G-W	O-NCG-QD-L-G-W	O-NCG500-QD-L-G	O-OX07-QD-L-XXX	O-OXSB-QD-L-XXX
DI	O-OH1-DI-L-G	O-OH2-DI-L-G	O-C73-DI-L-G-W	O-CSA-DI-L-G-W	O-CLI-DI-L-G-W	O-NCG-DI-L-G-W	O-NCG500-DI-L-G	O-OX07-DI-L-XXX	O-OXSB-DI-L-XXX
CH	Not Available	Not Available	Not Available	Not Available	Not Available	O-NCG-CH-L-G-W	Not Available	O-OX07-CH-L-XXX	Not Available
PB	Not Available	Not Available	O-C73-PB-L-G-W	O-CSA-PB-L-G-W	O-CLI-PB-L-G-W	O-NCG-PB-L-G-W	O-NCG500-PB-L-G	O-OX07-PB-L-XXX	O-OXSB-PB-L-XXX

**Using Existing Parts**  
 O-OH2-QD-L-G-PP  
 O-OH2-DI-L-G-PP

**L defines the Language:**

English (NFPA) = U  
 English (CSA/ISO) = E  
 French (CSA/ISO) = F  
 Spanish (NFPA) = S

PP refers to "plug and play" when using existing Diamond 2 & 3 rough-in

**G defines the Gas Type:**

Oxygen = O  
 Medical Air = A  
 Nitrous Oxide = 2  
 Carbon Dioxide\* = C  
 Nitrogen\*\* = N  
 MedVac = V  
 Waste Anesthetic Gas Disposal = W (NFPA)  
 Anesthetic Gas Scavenging System = AG (CSA/ISO)

\* Only Available in DISS or QD  
 \*\* Only Available in DISS

**W defines the Frame Width:**

3 = 3.25"  
 5 = 5.0"

**XXX defines the Gas Type:**

Oxygen = OXY  
 Medical Air = AIR  
 Nitrous Oxide = N2O  
 Carbon Dioxide\* = CO2  
 Nitrogen\*\* = NIT  
 MedVac = VAC  
 Waste Anesthetic Gas Disposal = WAG (NFPA)  
 Anesthetic Gas Scavenging System = AGS (CSA/ISO)

\* Only Available in DISS or QD  
 \*\* Only Available in DISS