


RIGID CEILING COLUMN

(M/N: C-COL-RIGD)


DRAWING # 60



TYPE:
QUANTITY:

SYSTEM DETAILS		
SYMBOL	QTY.	DESCRIPTION
CC	1	CEILING COLUMN
O	2	GAS, OXYGEN
A	2	GAS, MED AIR
V	3	GAS, VACUUM
2	2	GAS, NITROUS OXIDE
E	1	GAS, EVAC
E2	4	RECEPTACLE - DUPLEX RED
N2	2	RECEPTACLE - DUPLEX IVORY
TL	1	PROVISION - TELEPHONE 1G
SLM	1	VACUUM SLIDE MOUNTED

Two electrical/communication devices or gases cannot be placed together at any corner.



IMPORTANT: PLEASE VERIFY THAT THE ABOVE INFORMATION IS CORRECT, AND PROVIDE THE REMAINING DETAILS.

APPROVAL - PRINT AND SIGNATURE _____

DATE _____

PHONE NO. _____

AMICO Corporation
 85 Fulton Hwy,
 Richmond Hill, Ontario
 L4B 2N4, CANADA
 Tel: (905) 764-0980
 Fax: (905) 764-0982
 www.amico.com

HOSPITAL	HOSPITAL	A. NURSE CALL MFGR:	MODEL NO.:	DRWG. NO.
LOCATION	LOCATION	B. MEDICAL GAS MFGR.:	TYPE CONNECTION:	PRECBOOK-63
QTY.	1	C. FINISH:	CEILING HEIGHT:	DRAWN BY:
				CHECKED BY:
				REV. NO.:
				DATE: