


SELECT SERIES SECURITY WALL CONSOLE

DRAWING # 52


(M/N: W-X-SEC-PCON-07)
(M/N: E-CON-07)



TYPE: 7-GANG SECURITY CONSOLE
QUANTITY: 1

SYSTEM DETAILS	
SYMBOL	DESCRIPTION
SEC	SECURITY CONSOLE DOOR
O	AMICO GAS OUTLET - OXYGEN
A	AMICO GAS OUTLET - MED AIR
V	AMICO GAS OUTLET - VACUUM
SLM	VACUUM SLIDE - DOOR MOUNTED

- NOTE:
1. PUNCH OUT CENTRE TO CENTRE IS 2 1/2"
 2. THE BACK BOX WILL INCORPORATE 1/2" AND 3/4" KNOCKOUTS ON THE TOP AND BOTTOM. THE BACK BOX IS PROVIDED WITH ADJUSTABLE FLANGES TO ACCOMMODATE VARIOUS WALL THICKNESS.
 3. FASCIA WILL BE EXTRUDED ANODIZED ALUMINUM.



IMPORTANT: PLEASE VERIFY THAT THE ABOVE INFORMATION IS CORRECT, AND PROVIDE THE REMAINING DETAILS.

APPROVAL - PRINT AND SIGNATURE _____

DATE _____ PHONE NO. _____

AMICO Corporation
85 Fulton Hwy,
Richmond Hill, Ontario
L4B 2N4, CANADA
Tel: (905) 764-0800
Fax: (905) 764-0882
www.amico.com

HOSPITAL LOCATION: HOSPITAL LOCATION

QTY: 1 TYPE _____ UNITS AS SHOWN / (_____) TYPE _____ UNITS OPPOSITE

A. NURSE CALL MFGR.: _____ MODEL NO.: _____

B. MEDICAL GAS MFGR.: _____ TYPE CONNECTION: _____

C. FINISH: _____ CEILING HEIGHT: _____

DRWG. NO. _____ PRICEBOOK-89

DRAWN BY: _____ CHECKED BY: _____

REV. NO.: _____ DATE: _____