

Diagnostic Equipment

Aneroid Manometer



Speculum Holder



Sure Temp



Ophthalmoscopes



Features:

- More flexible positioning over the traditional wall board
- Equipment can be locked down or easily moved
- Compatible with all Welch Allyn 767 series
- Saves time on repairs, replacements and installation
- Can be mounted to any headwall, power column, flatwall, equipment boom or direct rail option

Innovative solutions that streamline workflow, improve patient safety and enhance patient care. Amico Accesories provides mounting solutions for all diagnostic equipment, assuring secure mounting and ergonomic solutions everywhere in the healthcare facility.



Mounting Options

- A** Amico Horizontal Rail Mount
(Compatible with Amico and Modular)
- V** Amico Vertical Rail Mount
(Compatible with Amico and Modular)
- F** Fairfield Rail Mount (push button optional)
- E** DIN 25mm Rail Mount (30mm, 35mm optional)
- C** 1 ¼" Pole Mount
- M** Monitor Rail Mount
- G** Vertical Rail Mount
(Compatible with Amico, Modular and Hill Rom)
- S** HSI Mount-Vertical Rail Mount

Amico Accessories solutions include the adapters required for the above rail systems.

Equipment Mounts

Amico BP Aneroid



- | | |
|--------------------|------------------|
| A AM-CB-A | G AM-CB-G |
| V AM-CB-V | S AM-CB-S |
| F AM-CB-F | |
| E AM-CB-E25 | |

*System includes; gauge, inflation system, adult cuff, and cuff basket

Cuff Basket



- | | |
|-------------------------|--------------------------|
| A AM-CUFBASK-A | C AM-CUFBASK-R125 |
| V AM-CUFBASK-V | M AM-CUFFBASK-MON |
| F AM-CUFBASK-F | G AM-CUFBASK-G |
| E AM-CUFBASK-E25 | S AM-CUFBASK-S |

*Mounting plate prevents cuffs from easily falling out

Welch Allyn 767 Aneroid



- | | |
|------------------------|-------------------------|
| A AM-767MNT-A | C AM-767MNT-R125 |
| V AM-767MNT-V | M AM-767MNT-MON |
| F AM-767MNT-F | G AM-767MNT-G |
| E AM-767MNT-E25 | S AM-767MNT-S |

*Mounting plate only, BP system not included

Welch Allyn 767 Ophthalmoscopes



- | | |
|----------------------------|-----------------------------|
| A OPTHALMNT-767-A | C OPTHALMNT-767-R125 |
| V OPTHALMNT-767-V | M OPTHALMNT-767-MON |
| F OPTHALMNT-767-F | G OPTHALMNT-767-G |
| E OPTHALMNT-767-E25 | S OPTHALMNT-767-S |

*Mounting plate only, Transformer not included

Welch Allyn 5091 Aneroid



- | | |
|-------------------------|--------------------------|
| A AM-5091MNT-A | C AM-5091MNT-R125 |
| V AM-5091MNT-V | M AM-5091MNT-MON |
| F AM-5091MNT-F | G AM-5091MNT-G |
| E AM-5091MNT-E25 | S AM-5091MNT-S |

*Mounting plate only, BP system not included

Welch Allyn Sure Temp or KleenSpec Holder



- | | |
|----------------------|-----------------------|
| A SPECHLD-A | C SPECHLD-R125 |
| V SPECHLD-V | M SPECHLD-MON |
| F SPECHLD-F | G SPECHLD-G |
| E SPECHLD-E25 | S SPECHLD-S |

*Mounting plate only, devices not included


Trimline 4300 Aneroid



- | | |
|-------------------------|--------------------------|
| A AM-4300MNT-A | C AM-4300MNT-R125 |
| V AM-4300MNT-V | M AM-4300MNT-MON |
| F AM-4300MNT-F | G AM-4300MNT-G |
| E AM-4300MNT-E25 | S AM-4300MNT-S |

*Mounting plate only, BP system not included


WA Baum Wall or ADC 750W Aneroid



- | | |
|-------------------------|--------------------------|
| A AM-750WMNT-A | C AM-750WMNT-R125 |
| V AM-750WMNT-V | M AM-750WMNT-MON |
| F AM-750WMNT-F | G AM-750WMNT-G |
| E AM-750WMNT-E25 | S AM-750WMNT-S |

*Mounting plate only, BP system not included

Please visit www.amico.com for all technical specifications, installation instructions, detailed AutoCAD drawings and more Instrument Mounting Solutions.



Distributed By:

