


Monitor roll stand



Features:

- Load capacity up to 25lbs (5.5kg)
- Overall height 42" (106.7cm)
- 3" (7.6cm) single wheel casters (2 locking, 3 non-locking)
- 5" (12.7cm) slide-in mounting plate or VESA 75/100 plate
- Custom color options (grey is standard)
- MRI safe options are available
- Sturdy 25" aluminum base

Part Number: RS1-LCD-42

Part Number: RS1-420-PKG

3" Casters

Patient Monitors	LCD Monitors
• Philips	• Toshiba
• Spacelabs	• Sony
• GE	• Samsung
• Nihon Kohden	• Gateway
• Datex Ohmeda	• Advantech
• Welch Allyn	• DELL
• Datascope	• NEC
• Zoll	• LG
Other Manufacturers Options Available	

Please specify model when ordering.

the right connection

www.amico.com

Roll Stand **Accessories**



Tablet / Laptop **Mount**

- Adjusts to fit most tablet/laptop PCs
- Tilt adjustment for easy viewing

Part Number: TABMNT-R125



Keyboard **Mount** (keyboard not included)

- Left or right mouse extension
- Tilt compatible for ergonomic typing in a sitting or standing position

Part Number: KEYBRDA-R125



Soft Grip **Handle**

- Easy grab handle

Part Number: RS1-HANDLE



Power **Bar**

- Hospital grade, 4 power outlets
- 120 VAC, 60Hz

Part Number: POWER4-R125



Power **Holder**

- Holds most power adapters for PCs, patient monitors and LCDs

Part Number: PWRHLD-R125



Cord **Wrap**

- Conveniently wrap cords and keep them organized

Part Number: CORDWRAP-R125



Wire **Basket**

- 10"W x 6"D x 9" H white wire basket
- Convenient storage

Part Number: 100609-R125-BB



Cable **Manager**


- Easy pass through channel to manage multiple cables

Part Number: MRS-009-R125

Please visit www.amico.com for more information on our complete line of mobile solutions.



Distributed By:



SP-RS1062

the right **connection**