

## SSM monitor arm

SIDE TO SIDE MOTION



### Features:

- 9" or 12" (23cm or 30cm) single arm lengths
- Available for wall channel or rail mounting
- Swivel arm provides lateral functionality
- Tilt and swivel adjustments
- 60 lbs (27.2kg) arm rating
- Integrated cable management
- Color options available (white is standard)



Fully adjustable tilt and pivot

## Monitor Compatibility

Patient Monitors	LCD Monitors
• <b>Draeger</b>	• <b>Toshiba</b>
• <b>Spacelabs</b>	• <b>Sony</b>
• <b>GE</b>	• <b>Samsung</b>
• <b>Nihon Kohden</b>	• <b>Gateway</b>
• <b>Datex Ohmeda</b>	• <b>Advantech</b>
• <b>Welch Allyn</b>	• <b>DELL</b>
• <b>Datascope</b>	• <b>NEC</b>
• <b>Philips</b>	• <b>LG</b>
<b>Other Manufacturers Options Available</b>	

Please specify model when ordering.



Part Number: SSM-LCD-15  
12" x 12" Dual Arm

the right connection

www.amico.com

# Additional **SSM Models**

## LCD Flush Mount



- 2.25" (6cm) profile
- 30 lbs (13.60kg) arm rating
- VESA 75/100 plates included
- Portrait and landscape viewing


Part Number: SSM-LCD-45

## Flush Mount



- 2.25" (6cm) profile
- 30 lbs (13.60kg) arm rating
- Back mount devices only

Part Number: SSM-456



## LCD Double **Pivot Arm** 9"x 9"

- 16" (41cm) extension
- 30 lbs (13.60kg) arm rating
- VESA 75/100 plates included
- Portrait and landscape viewing

Part Number: SSM-LCD-35



## Double **Pivot Arm** 9"x 9"

- 16" (41cm) extension
- 30 lbs (13.60kg) arm rating
- ± 30° degrees tilt

Part Number: SSM-350



## LCD Double **Pivot Arm** 12"x 12"

- 22" (56cm) extension
- 30 lbs (13.60kg) arm rating
- VESA 75/100 plates included
- Portrait and landscape viewing

Part Number: SSM-LCD-15



## Double **Pivot Arm** 12"x 12"


- 22" (56cm) extension
- 30 lbs (13.60kg) arm rating
- ± 30° degrees tilt

Part Number: SSM-150

Please visit [www.amico.com](http://www.amico.com) for more information on our complete line of monitor accessories.



## Distributed By:



SP-SSM062

the right **connection**