

INNOVATION THAT MAKES SENSE

Integrated Flowmeter...

Gas Specific Labeling and Fittings
Medical Oxygen and Medical Air

Supply Outlet 50 PSI

Adjustable Knob for flow control
Ranges available:
0-8 LPM, 0-15 LPM, 0-25 LPM

Flow Indicator

Flowmeter Port

- ✓ 2 outlets in 1
- ✓ Easy to service
- ✓ Retrofit existing outlets
- ✓ No more lost, chained, or broken Flowmeters
- ✓ Safest, most convenient way to dispense Medical Oxygen & Air

SP **RU**

Amico Integrated Flowmeter outlet combines a Dial Flowmeter and a Medical Gas Outlet into a single compact design. The right side port connection has a control knob which allows the user to adjust the flow setting, while the left side port is a direct connection to the gas supply.

The Amico Integrated Flowmeter saves space, reduces construction costs, and ensures healthcare providers always have a Flowmeter when they need one.

model numbers

S-FI-CXXYZ-WQ-S

Flow Range:

8 LPM: 08
 15 LPM: 15
 25 LPM: 25

Colour Code:

USA: U
 ISO: I

Gas Type:

Oxygen: O
 Medical Air: A

Line Fitting:

DISS Male with Check Valve: D
 Ohmeda Coupler: O
 Chemetron Coupler: C

Flowmeter Fitting:

DISS Male with Check Valve: D

Amico Retro-Fit Conversion Compatibility:

Diamond 1 & Diamond 2/3 (Complete Conversion):	A
Diamond 2/3 (Plug and Play):	B
Puritan Bennett Series 73, Series A & Litharge NCG/Chemetron Series 500, 560	
Oxequip/Medstar Series 27 (Smooth Bore):	C
NCG/Chemetron Series 148, 400, 460 Oxequip/Medstar Series 7:	Leave Blank
(Regular) Amico Rough-in:	Leave Blank

NOTE: Rough-in conversions, Filler Frames and Trim plates are not included in these part numbers.

represented by: