

Our Medical Gas Outlets are designed to update and convert your old Outlets.

The Retro-Fit process allows you to update and convert your old Medical Gas Outlets to a new Amico Outlet. This can be done without cutting the pipeline or breaking into the existing wall. Amico's "Retro-Fits" are designed to use the existing back body of the previous product. This allows you to maintain the original configuration, while improving both the functionality and aesthetic appearance of the Outlet.

Old

Diamond 1 Outlet



New

Amico Ohmeda Compatible or DISS Outlet



Old

Daimond 2 Outlet



New

Amico Ohmeda Compatible or Diss Outlet



Old

Daimond 3 Outlet



New

Amico DISS or Ohmeda Compatible Outlet



project

outlets - retro-fit

Old



New




Old



New



Old



New



represented by: